



UPGRADE YOUR AIRWAY SIMPLIFY YOUR WORKFLOW

For too long, anesthesia providers have improvised ways to achieve adequate oxygenation, capnography, and airway support during deep IV sedation. AERIS™ provides an elegant, simple, FDA-cleared solution by converting standard OPAs and NPAs into an **Advanced Sedation Airway Platform**—delivering high-efficiency oxygenation and dependable capnography to bridge the gap before airway escalation.

See How It Works:



Designed for
Deep IV Sedation



Patients Who Benefit the Most

AERIS™ is designed for challenging deep IV sedation patients, including those with:

- ▶ Higher BMI
- ▶ OSA / mouth breathing
- ▶ Challenging airway anatomy
- ▶ Restricted access to the face
- ▶ Longer cases




Add the AERIS Airway™ to your OPA or NPA

to provide the support needed to avoid airway interventions and more invasive airway management.

The Airway Continuum: Support at Every Step

MINIMAL SUPPORT Nasal Cannula / Face Mask



Limited oxygen delivery
and no airway stabilization.

ADVANCED SUPPORT




+ standard
OPA / NPA



- Enhanced oxygenation
(up to 80-90% FiO₂ at 2-4L/min).
- Real-time, reliable capnography.
- Supports airway patency.


BRIDGES THE GAP BEFORE
AIRWAY ESCALATION.

MORE INVASIVE LMA



Requires deeper sedation
or general anesthesia.
Escalated airway management.

MOST INVASIVE ETT



Requires general anesthesia
with mechanical ventilation.



"AERIS™ provides much more effective oxygenation than nasal cannulas and delivers an end-tidal CO₂ waveform during spontaneous ventilation that is remarkably close to that of an endotracheal tube. This is an amazing device."

ANDREW HERLICH, DMD, MD, FAAP, FASA
DUAL-TRAINED IN DENTAL AND MEDICAL ANESTHESIOLOGY; PROF. EMERITUS AND CHAIR OF THE DEPARTMENT OF ANESTHESIOLOGY, UNIV. OF PITTSBURGH



"The AERIS Airway™ takes the guessing and estimating out of my practice when monitoring an anesthetized patient. AERIS™ provides confidence that CO₂ monitoring is within ASA guidelines."

DIANA LUGO, DNAP, CRNA
FORMER PRESIDENT, ARIZONA ASSOCIATION OF NURSE ANESTHESIOLOGY




"AERIS™ administers O₂ and monitors CO₂ significantly better than divided nasal cannulas and face masks."

RENE M. GONZALEZ, MD
FORMER ASSOCIATE PROF. OF ANESTHESIOLOGY AND CRITICAL CARE MEDICINE, UNIV. OF PITTSBURGH SCHOOL OF MEDICINE

AERIS™ Advantage


Transforms Oral and Nasopharyngeal Airways into State-of-the-Art O₂ Delivery & CO₂ Monitoring Systems



AERIS Airway™ is compatible with:
OPAs (80-110mm)
NPAs (24-36Fr)


O₂ Delivery Port:

- O₂ delivery to pharyngeal region
- Superior FiO₂ vs. Nasal Cannula and Face Mask; approaching 80% at flow rates as low as 2-4 L/min
- Saves on O₂ costs and minimizes fire risk



CO₂ Detection Port:

- Captures real-time, reliable CO₂ monitoring
- Enable anesthesia provider to preemptively detect & prevent respiratory compromise
- Reduces the risk of hypoxemia



Ref #	Description	UOM
1030MM	AERIS Airway™ Male-Male Luer Kit*	30/Bx
1030MF	AERIS Airway™ Male-Female Luer Kit*	30/Bx
1030MSA	AERIS Airway™ Male-MSA Luer Kit*	30/Bx

*AERIS Airway Kit contains: AERIS Airway™, 0ft Nasal Cannula, 10ft bonded O₂/CO₂ tubing