







3 AERIS™ Kits Designed to Fit All Existing Patient Monitors

Ref#	Description	Qty/Box	Patient Monitor CO ₂ Luer Connection**
1030MM	AERIS Airway™ Male-Male Luer Kit* 	30/box	Male 
1030MF	AERIS Airway™ Male-Female Luer Kit* 	30/box	Female 
1030MSA	AERIS Airway™ Male-MSA Luer Kit* 	30/box	MSA 

* AERIS™ Kit contains: AERIS Airway™, 0ft Nasal Cannula, 10ft bonded O₂/CO₂ tube set

** Existing Patient Monitors use either a Male, Female, or MSA CO₂ Luer connection