

Introducing

# AERIS AIRWAY™

THE FUTURE OF AIRWAY MANAGEMENT



AERIS™ in OPA



AERIS™ in NPA



## SIMPLE.

The AERIS Airway™ seamlessly integrates across existing oropharyngeal airways (OPA) and nasopharyngeal airways (NPA).

## SAFER.

AERIS Airway™ delivers maximal FIO<sub>2</sub> and a real-time, reliable CO<sub>2</sub> waveform, allowing anesthesia providers to preemptively detect and prevent hypoxemia.

## SEDATION.

A first-of-its-kind oxygen delivery and end-tidal CO<sub>2</sub> monitoring device for I.V. sedation.



*"The AERIS Airway™ administers O<sub>2</sub> and monitors CO<sub>2</sub> significantly better than divided nasal cannulas and face masks."*

**RENE M. GONZALEZ, MD**  
FORMER ASSOCIATE PROF. OF ANESTHESIOLOGY AND CRITICAL CARE MEDICINE, UNIVERSITY OF PITTSBURGH SCHOOL OF MEDICINE



*"AERIS™ is a game changer for patient safety and monitoring ETCO<sub>2</sub>."*

**ERIC WARDRIP, MD**  
MEDICAL DIRECTOR, SCA CENTER FOR SURGERY OF ENCINITAS




*"The AERIS Airway™ is a breakthrough in airway management, enabling anesthesiologists to preemptively detect signs of compromised airway flow and prevent hypoxemia."*

**ALTON LIU, MD**  
PARKVIEW REGIONAL MEDICAL CENTER


## AERIS™ Advantage

Transforms Oral and Nasopharyngeal Airways into State-of-the-Art O<sub>2</sub> Delivery & CO<sub>2</sub> Monitoring Systems




### O<sub>2</sub> Delivery Port:

- O<sub>2</sub> delivery to pharyngeal region
- Superior FIO<sub>2</sub> vs. Nasal Cannula and Face Mask; approaching 80% at flow rates as low as 2-4 L/min
- Saves on O<sub>2</sub> costs and minimizes fire risk



### CO<sub>2</sub> Detection Port:

- Captures real-time, reliable CO<sub>2</sub> monitoring
- Enable anesthesia provider to preemptively detect & prevent respiratory compromise
- Reduces the risk of hypoxemia



AERIS Airway™ is compatible with:  
OPAs (80-110mm)  
NPAs (24-36Fr)



[Click here to see how the AERIS Airway™ works](#)

Ref #	Description	UOM
1030MM	AERIS Airway™ Male-Male Luer Kit*	30/Bx
1030MF	AERIS Airway™ Male-Female Luer Kit*	30/Bx
1030MSA	AERIS Airway™ Male-MSA Luer Kit*	30/Bx

\*AERIS Airway Kit contains: AERIS Airway™, 0ft Nasal Cannula, 10ft bonded O<sub>2</sub>/CO<sub>2</sub> tubing